

Nextractor® NX-48S Platform for COVID-19 Viral RNA Extraction

Viral RNA Extraction in 17 min for COVID-19 Detection

| Feature | Benefit |
|--------------------------------------|---|
| High throughput | 48 specimens per run |
| Fast extraction time (17 min) | Quick COVID-19 diagnosis |
| Fully automated extraction | Minimizes hands-on time |
| Compact size | Fits in BSC Multiple unit installation possible in small space |
| UV sterilization | Eliminates contamination |

Extraction Instrument

 **Nextractor® NX-48S**



Extraction Kit

| Product | Package Size | Cat. No | Specimen Type |
|--------------|-------------------|---------------|---------------------------------------|
| Viral NA Kit | 96 test / 48 test | VN142/152/162 | Raw sputum |
| | | VN143/153/163 | Nasopharyngeal swab Diluted sputum |